

Setting Up the 2D Spectral Window

- Resolution is Limited in Both Dimensions, but Especially in F_1 (Indirect) Dimension
 - 1D FID has 16,384 or 32,768 Data Points
 - F_2 FID of 2D Spectrum has 512 to 2048 Points
 - F_1 FID of 2D Spectrum has 128 to 750 Points
- “Wasted Space” in Spectral Window Must Be Eliminated
 - Include Only Peaks from Your Compound
 - ^1H Spectrum: Cut Out TMS and CHCl_3
 - ^{13}C Spectrum: Cut Out Quaternary Carbons for HSQC
- Spectral Window is Narrowed in F_1 and F_2
- Bruker Parameters For Setting Up Spectral Window:
 - **sw**(F_2), **sw**(F_1): Spectral Width in ppm
 - **o1p**: Center of Window in ppm (F_2)
 - **o2p** or **o1p**: Center of Window in ppm (F_1)
 - For Heteronuclear 2D: F_1 corresponds to **o2p**
- Example: ^1H : 1.0 to 6.0 ppm; ^{13}C : 20 to 140 ppm
sw(F_2): 5.0 ppm **o1p**: 3.5 ppm
sw(F_1): 120 ppm **o2p**: 80 ppm